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**realme**  
realme Smart Scale  
Operating Manual  
RMH2011

**Warning**

- The product is forbidden to be dismantled, extruded, stabbed or short-circuited randomly or put in water or fire or exposed in an environment with a temperature higher than 40°C (104°F).
- A person installed with a cardiac pacemaker or other medical appliances is forbidden to use this product. A pregnant woman should use this product with caution and should follow the doctor's advice.
- CAUTION:** Risk of explosion if the battery is replaced by an incorrect type. Replace only with the same type alkaline battery.
- Dispose of a battery into fire or a hot oven, or mechanically crushing or cutting of a battery, that can result in an explosion.
- Leaving a battery in an extremely high temperature surrounding environment that can result in an explosion or the leakage of flammable liquid or gas.
- A battery subjected to extremely low air pressure that may result in an explosion or the leakage of flammable liquid or gas.

**Important Instruction**

- This product is not a medical appliance.
- The measured data from the realme body fat scale are for reference only and should not serve as a basis for medical opinions.

**Declaration**

Body data about the following groups may be inaccurate:

- Minor
- Professional athlete
- Pregnant woman
- A patient with a cardiac pacemaker

**Basic Parameters**

- Product Model: RMH2011
- Product Size: 300\*300\*23.25mm (including the electrode and scale base)
- Net Weight of the Product: 1.74kg
- Unit: kg/lb
- Display Method: LED panel
- Power Source: 4\*1.5V AAA batteries
- Storage environment: -20°C ~ 60°C
- Working environment: 5°C ~ 40°C
- Frequency Range: 2402 ~ 2480MHz
- Transmit Power (EIRP): 0.2dBm

**Packing List**

- realme Smart Scale: 1
- 1.5V AAA battery: 4
- Operating manual & warranty card: 1

**Product Information**

**Operating Manual**

**1 Install batteries**

Take the scale from the plastic bag, open the battery case at the bottom of the scale, pull out the insulation sheet from the battery jar and cover the battery cover.

**3 Bind the body fat scale**

Open the realme Link APP and click "Go to the top right corner of the home page."

Enter the page of realme Smart Scale and bind it.

Choose realme Smart Scale and bind it.

**4 Click the button of "Step on" in the APP and measure yourself barefooted**

Begin to use

- Correct standing posture
- The feet should be distributed naturally at both sides with the scale center as an original point

**5 Heart rate measurement: Heart rate measurement will begin after the data about weight is locked and displayed for 3 seconds.**

When the "Step on" button in the top right corner of the screen, the body fat scale will enter the mode of heart rate measurement.

When the measurement is going on, an animation of feet bounding of people in the middle of the screen.

After the measurement is completed, the animation will stop and the data about heart rate will be locked and appear on the screen.

**Troubleshooting**

- The mobile phone configuration does not support the product. Please refer to the following requirements about mobile phone system Supporting system: Android 5.0 and above, Bluetooth 4.0 and above
- Insufficient electric quantity of the scale
- After the battery voltage is lower than 4.2V, the body fat scale will initiate lower electricity reminder and the screen will show "Lo". But the use of the product will not be affected. When the battery voltage is lower than 4.0V, the body fat scale will enter a state of empty power and the screen will show "Lo". The equipment will be automatically powered off in 3s. The screen will show "Lo". The equipment will be automatically powered off in 3s. The screen will show "Lo". The equipment will be automatically powered off in 3s. The screen will show "Lo".
- The network signal for the mobile phone is weak
- (WiFi or 4G network is recommended)
- The mobile phone Bluetooth is not initiated

**Warning**

**1. EU Declaration of Conformity (DoC)**

Hereby, realme, declares that this device is in compliance with the essential requirements and other relevant provisions of Radio Equipment Directive 2014/53/EU. The declaration of conformity may be consulted at [www.realme.com/global/support/eu-declaration](http://www.realme.com/global/support/eu-declaration).

**2. Recycling Initiative**

At realme we understand that our responsibility doesn't end at selling you our products. realme have been working in the area of safe disposal of electronic wastes. Since mobiles, wearables and other electronic items are made of hazardous constituents. Under the Guidelines of Ministry of Environment, Forest and Climate Change, Government of India, E-waste (Management) Rule, 2016 and Amendment E-waste Rule 2018:

- realme will seek shared responsibility and cooperation from customers in reducing the environmental impact of their products.
- realme will comply with all the applicable laws related to e-waste management.

**Warranty Policy and Customer Service Information**

Thank you for purchasing the product of realme. If any quality problems arise within 12 months from the date of purchase, you can enjoy the company's warranty service.

**Matters Needing Attention:**

- When you buy the product, the sales unit will issue a valid proof of purchase.
- The warranty card and proof of purchase should be kept properly and will not be reissued if you lose them.
- As to failures caused by non-human factors, the user can enjoy free maintenance at the customer service center of realme by showing his/her warranty card and proof of purchase within the warranty period.

**Non-Warranty Circumstances:**

Any of the following cases is not within the warranty scope during the using process.

- The warranty period expired
- Damages caused by using the product not in accordance with the operating manual or operating guide.
- Damages caused by force majeure factors (flood, earthquake, lightning, etc.)
- There is no warranty card or invoice or the information of the warranty card is inconsistent with the information of the invoice.
- Natural abrasion of the product.
- Failures and damages caused by factors other than quality problems of the product.

**Warranty Card (Stub Copy)**

Dear User, thank you for using the product of this company. Please read and keep properly the warranty card after purchasing the product so that we can provide better services to you.

**User Information**

User's Name \_\_\_\_\_ Purchase Date \_\_\_\_\_  
 Contact Tel \_\_\_\_\_ Invoice Number \_\_\_\_\_  
 User's Address \_\_\_\_\_ Sales Unit \_\_\_\_\_  
 Email \_\_\_\_\_ Sales Address \_\_\_\_\_

**Product Information**

Product Model \_\_\_\_\_ Tel \_\_\_\_\_  
 Product Serial Number \_\_\_\_\_

**Warranty Card (Customer Copy)**

Dear User, thank you for using the product of this company. Please read and keep properly the warranty card after purchasing the product so that we can provide better services to you.

**User Information**

User's Name \_\_\_\_\_ Purchase Date \_\_\_\_\_  
 Contact Tel \_\_\_\_\_ Invoice Number \_\_\_\_\_  
 User's Address \_\_\_\_\_ Sales Unit \_\_\_\_\_  
 Email \_\_\_\_\_ Sales Address \_\_\_\_\_

**Product Information**

Product Model \_\_\_\_\_ Tel \_\_\_\_\_  
 Product Serial Number \_\_\_\_\_

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